

Supplementary Worksheet 1 - BioHansel

Geneviève Labbé et al

2019-11-26

Detailed parameter settings for all programs used for validation

Concordance Assessment - *S. typhi*

Command-lines for BioHansel and Snippy found in Nextflow workflow nf-biohansel-snippy-comparison

BioHansel command-line:

```
mkdir -p reads
ln -s `realpath *.fastq.gz` reads/
hansel \\\
-v \\\
-t 1 \\\
--min-kmer-freq $mincov \\\
-s $scheme_fasta \\\
-D reads/ \\\
-o $summary_report \\\
-O $detailed_report
```

Parameters:

- -v verbose logging
- -t 1 run BioHansel with one CPU thread
- -D reads/ run BioHansel on reads symlinked to reads/ directory
- \$scheme_fasta path to *S. typhi* BioHansel scheme
- \$mincov minimum kmer frequency value of 3, 6, 8
- \$summary_report BioHansel summary report output file
- \$detailed_report BioHansel detailed results output file

Snippy command-line:

```
snippy --prefix $accession \\\
--outdir $accession \\\
```

```
--cpus 16 \\  
--ram 16 \\  
--mincov $mincov \\  
--R1 $reads1 \\  
--R2 $reads2 \\  
--ref $ref_genbank \\  
--tmpdir ./
```

- \$accession SRA run accession (sample name)
- \$mincov minimum coverage depth of 3, 6, 8
- \$reads1/\$reads2 forward and reverse paired-end Illumina reads paths
- \$ref_genbank GenBank file for reference genome NC_003198.1 (*Salmonella enterica* subsp. *enterica* serovar Typhi str. CT18, complete genome.)